

New upgrade in 2020

808nm laser hair removal upgrade

- Full absorb
- 808+532+755+1064nm
- Suitable for Asians
- Non-invasive and painless

For hair removal machines, hair removal is more professional and safer.



Are you often troubled by hair?



Not daring to lift an arm



Dare not wear short skirts



Afraid to wear a bikini



8 weeks smooth

Natural shedding / painless hair removal / no black spots



Delicate
skin



8
Week



1-4 weeks

2-3 times a week, the hair began to fall off naturally, and the hair was reduced by 75%.

5-6 weeks

2 times a week, new growth hair falls off and inhibits growth, maintaining a smooth skin.

7-8 weeks

Once a week, the hair becomes thinner, less, no longer grows, and keeps the skin smooth.

More than 8 weeks

1-2 times a week, depending on the state of hair growth, selective use, you can say goodbye to repeated hair removal.

Spot size:15mm*15mm



SCOPE OF APPLICATION

////////// Efficient //////////



Armpit hair removal



Facial hair removal



Beard removal



Chest hair removal



Arm hair removal



Chin & lip hair removal



leg hair removal



Bikini hair removal

Before and after use

Healthy hair removal, the effect is visible



BACK HAIR REMOVAL



LEG HAIR REMOVAL



ARMPIT HAIR REMOVAL



CHIN & LIP HAIR REMOVAL



Before

After



Before

After



Before

After



Before

After

LASER BEAUTY HANDLE

////////// New two-in-one //////////



Holographic focusing technology can focus the high-energy laser on the skin to stimulate collagen proliferation. Go blackheads, wash eyebrows, tattoos, freckles.

Black Baby tip

Removes fine wrinkles, dilutes true wrinkles, tightens skin, shrinks pores, removes blackheads and acne, dilutes acne marks, controls oil, and improves oil secretion. It can prevent acne and treat early acne.





1064 gun head

Mainly black, blue, purple with darker pigments. (Birthmarks, Dermal Spots, Moles, Age Spots, Ota Moles, Ito Moles, Tattoos, Eyebrows, Eyelines.)



532 gun head

Lorem ipsum dolor sit amet, consectetur adipiscing elit, sed do eiusmod tempor incididunt ut labore et dolore magna aliqua. Quis ipsum suspendisse ultrices gravida. Risus commodo viverra maecenas accumsan lacus vel facilisis.



Black Baby tip



Freckle



Wash eyebrows



Tattoo removal



PRODUCT DETAILS

////////// Everywhere details //////////



• *Foldable operation screen*





Diode laser parameters

Wavelength	808nm
Channel	macro channel
Power supply	2000W
Output power	300W (600W optional)
Energy	1-120mj
Frequency	1-10Hz
Pulse width	1-400ms
Cooling	air+water+semi conductor

Nd yag laser parameters

Wavelength	532nm 1064nm 1320nm
Energy	1-2000mj
Frequency	1-10hz
Power supply	800W

Machine size	44cm*44cm*133cm (the display stands up) 44cm*44cm*109cm (operation screen down)
--------------	--

Screen display





