

# 5 in 1 Ultrasonic 40K Cavitation & RF Slimming Unit



### 40K Cavitation Head

Power supply output : 150v  
Frequency : 400 Watt  
Power : 25 Watt



### Vacuum + Bipolar RF Head

Power supply output : 66v  
Frequency : 1 MHz  
Power : 95 Watt



### Tripolar RF Head

Power supply output : 66v  
Frequency : 450 khz  
Power : 35 Watt



### Face RF Head

Power supply output : 66v  
Frequency : 450 khz  
Power : 35 Watt



### Body RF Head

Power supply output : 66v  
Frequency : 550 khz  
Power : 25 Watt



<b>Energy output</b>	56 x 180 mW
<b>Main power output</b>	100 VA
<b>Safety</b>	On board diagnostics
<b>Classification</b>	Class 1, Type B applied part
<b>Cooling requirements</b>	Air cooled
<b>Mode of operation</b>	Continuous



## Functions:

Firming the arms, legs, thighs, buttocks, waist, back, abdominal muscles and reshape the body figure.

Improve and treat the orange peel like fatty tissues of different serious levels with the effects of soothing, decompression and pain relief (such as joint pain, sciatica), promote the regular mechanical aerobic exercise of skin tissue, rapidly relieve stress and remove fatigue.

Firming, lifting and compact the skin and mould the face shape circulate and activate cells, promote metabolism and so that 100% more collagen and elastic fibers may be gained for improving the facial and body skin.

## Features

1. Painless treatment focus RF energy at the correct position, compared to other radio frequency technology. It uses low energy and high frequency, safe and effective.
2. Aim at the skin surface and deep location, using a complex method to control the different currents and energy, direct access to the different skin layer. Does not cause the phenomenon of uneven skin.
3. Selectively targeting fatty tissue, to avoid heating other fat, to achieve the fastest treatment effect.
4. The entire process is completed without the need of surgery and anesthesia.
5. With the most useful 40KHZ ultrasonic cavitation system.
6. Without side effects and the risk of gaining the weight back, the results are rather obvious. Does not affect the normal work and life.

## Handles included:

### 40K Cavitation Head:

Cavitation causes vibration of the liquid parts inside a fat cell, the vibration erupts the fat cell from the inside out, exposing the fat, and eliminating the fat pocket and membrane.

### RF (Radio Frequency) Handle:

Radio frequency tones, tightens and stimulates collagen in the skin. Lifting and firming of areas such as the buttocks, thighs, breasts, face and eyebags. This unit has a small, medium and large handle.

### Vacuum Handle:

Lifts areas and is used to soften the appearance of cellulite.

**R15 999.00**

**Slimming technology that works**