

CE & EC Certification product



Recovery Technology Physiotherapy Technology

The best professional choice for private or commercial use

- ✓ Muscle Pathology
- ✓ Bone Injuries
- ✓ Inflammation and nerve damage

Terapia Magnetica Dual Pulsed System

Used mostly in physiotherapy healing, magneto Therapy pulses a magnetic field into the body, creating an extraordinary healing effect. The results are less pain, reduction in swelling, and increased range of motion in the affected areas. Each session lasts 5 -20 minutes, and a treatment cycle is 6 to 8 sessions depending on severity.



Specifications:

Power	1000 watt and 220 Volts
Pulse Continuous Time:	≥250 us (±20us)
Magnetic induction	2 Tesla
Laser	50watt
Wavelength	685nm- 830nm

System weight	60 kg
	57cm Length x 49cm Width x 113cm Height
Applicators	Double handle
Movement	4 lockable castor wheels