
3-in-1 Next Gen SLIMMING System 635nm Red light & 532nm Green light







635nm Red Light



532nm Green Light





Cryolipolisis is an unique non-invasive method for the selective reduction of fat cells with controlled localized cooling. This fat freezing Crystallization procedure makes fat cells under a process "apoptosis". Application delivers a great comfortable cooling treatment without oedemas.

How to use:

Step 1:

10 laser beams with continues or pulsed setting. 635nm Red Light Laser and 532nm Green Light laser.

Step 2:

4 cryo pads, with tempratures range between -10 to 4 degrees celcius.

Step 3:

Switch the pads to Electro Muscle Stimulation mode, to tone and tighten the muscle in the area


Why Choose The 3-in-1 Next Gen Slimming System?

- l. Market Approved to Overall Body Circumference Reduction
- 2. Highest success rate of all other laser fat treatment
- 3. Amazing results in as little a 2 weeks
- 4. No set-up time, can be used unattended as it is not applied directly onto the skin
- 5. Irradiation by a specific wavelength of cold source laser
- 5. Shrinks fat cells instead of killing them making it safer & more effective
- 7. No bruising, no freezing, no swelling, no downtime.

Specifications

Nett Weight Gross Weight	39 kg 58 kg	Screen size	12.1 Inch Touch Screen
Machine Size	118x58x60cm	Cooling System	Air Cooling
Treatment Area	30x800mm	Intensity	1 - 3 Level
Light Panel	10 Pcs		•
Power	800 Watt		

R195 000.00



Rejuvenate and refresh your appearance by tightening, lifting, stimulating collagen and refining skin texture for a smooth finish