## VShape<sup>3</sup> Cellulite and conturing system

VShape<sup>3</sup> is the gold standard for the treatment of cellulite and circumferential reduction

80% Of all women over the age of 20 have cellulite, regardless of size or shape. Cellulite occurs when fat cells in the deepest layer of the skin become over-inflated from an increase in fat storage.

The fat cells protrude up through the tissue causing a dimpling or a "cottage cheese" appearance. It uses a combination of bipolar radio-frequency, infra-red, and massage to achieve the desired results.


A course of 6 to 8 sessions once a week will produce the best results and each session takes on average 40 minutes and you can return to daily life straight afterward.



Infrared RF Vaccum Roller Slimming Machine is a super device, which gathers the technology of vaccum negative pressure, mechanicle rolling and skin folding, applying to all parts of shaping and wrinkle reduction.

## **Features**

2 Vacuum handsets
Precise massage roller technology
Infrared laser technology
Ultrasound cavitation
Bi-polar RF technology



Screen	10.4 inch colour touch screen	
Handle number	4 heads available	
Handle screen	2.8 inch & 3.5 inch	
Handle action area	4.5mm x 9mm x 13mm	
	16mm x 25mm (small handle)	
	30mm x 44mm (mid handle)	
	40mm x 120mm (big handle)	
	90mm x 120mm (cavitation har	ıdle)
Vacuum	10kpa - 75kpa	
Pulse width	1s - 8s	
IR wavelength	940nm	
Ultrasonic frequency	40kHz	
RF frequency	1Mhz	
Power supply	AC220V, 50Hz/AC110V, 60Hz	
1	I .	

A non-invasive treatment for face and body contouring, and master choice for cellulite treatment