

10-in-1 HOCATT™ PLUS

When it comes to the HOCATT™ Plus Ozone Technology, we in the wellness industry have a saying: Built for Wellness. Backed by Science. This is more than a bold claim. It's a statement of fact. The HOCATT™ Plus has 10 different technologies that were developed with a single purpose in mind – strategic wellness and vitality. It provides some of the top alternative health and beauty modalities in the world! That's why the HOCATT™ Plus is one of the most widely used systems globally.

**The ultimate wellness system
that delivers results that
matter.**



R244 662.50

Concentrates oxygen from the air in 93 - 97% pure oxygen. Feeds concentrated oxygen to the HOCATT's two ozone generators and also passes oxygen through the humidified oxygen breathing apparatus, so concentrated humidified oxygen may be inhaled during the HOCATT session.



Technical Specifications

Chamber

Voltage	220-240VAC, 7A, 50Hz1000W
Power	1500W
Weight	110kg
Dimensions	140x130x67cm
Heart Rate Monitor	Yes

Medical Concentrator

Voltage	220-240VAC, 2.2A, 50Hz
Power	480 W
Dimensions	6.4x12.7x30cm
Flow Rate	1 - 5L /m

Additional Features

Luxury Smith Chamber (fiberglass chamber with special ABS coating)
Auxiliary Ozone Functions can also be used separately
1 x Adjustable Seat
1 x Adjustable Foot Rest
1 x No1 Ozone Generator (for ozone sauna)
1 x No2 Ozone Generator (for simultaneous auxiliary ozone functions)
1 x FSM Foot Plate Set and 1 x FSM Electrode Set
1 x Electric CO ² Regulator
1 x Water Resistant Heart-Rate Monitor
1 x Full Colour Touch Screen Controller
1 x Reusable Oxygen Cannula
1 x Oxygen Humidifier
4 x Flow Meters

Accessories and Consumables:

5 x Neck Capes
3 x Glass Funnels (for cupping)
3 x Breast Cupping Set (S,M,L)
1 x Bag of Disposable Vaginal Catheters (50pcs)
1 x HOCATT Silicone Tube Set
3 x Water Dispensing Bottles

Features

Transdermal Ozone

When Ozone is introduced into the blood, it reacts with oxidizable substrates, killing viruses, bacteria and fungi, enhances blood circulation, and stimulates the Oxygen metabolism of the body.

Micro Currents

Micro Currents is an FDA approved Technology that can be applied in a variety of treatment modalities, Micro currents can assist with a wide variety of conditions, including pain and aches, nervous system, mental and physical, wellbeing & more

Insufflation

VI is recommended for endometriosis, ovarian cysts and vaginal infection (yeast/fungal, bacterial, viral). Often these infections decrease fertility, so VI is highly recommended for those who are looking to get pregnant.

Far Infrared

Far Infrared is a low energy Technology where waves of light penetrate the human body, transforming into heat energy to gently elevate temperature and activate major bodily functions. Benefits include: reduced muscle & joint pain, lower hypertension, sooth digestive issues, & much more.

Hyperthermia

Hyperthermia - mimics a fever state resulting in an effective natural process of detoxifying and stimulating healing processes, with all the same healing effects as the natural fever.

Photon Color Light

This state-of-the-art digital lighting Technology (LEDs) generate a variety of coloured lights to which the body reacts. In the HOCATT two colour stars change randomly to display multicolour cross fades; from soothing and relaxing to balancing & energizing.

Bagging

A limb is covered with an airtight bag and concentrated ozone is introduced to treat the wounds.

Essential Oil Infusions and Aromatherapy

The Essential oils are infused with the steam as it fills the chamber, making this a transdermal application. Used together has many health benefits.

Oxygen Respiration

Pure, humidified, oxygen runs through a special oxygen jet to the nose and mouth, boosting Oxygen saturation in the body during the HOCATT™ session.

CO2 I Carbonic Acid

Carbonic Acid is better known for its potent vasodilation ability of the skin blood vessels, which stimulate circulation, and oxygenation of the skin. For this reason, this modality is popularly used in cosmetic and reconstructive surgery centres and high-end spa's and beauty salons.

Cupping

Using a glass funnel/set of breast cups focuses a high concentration of ozone over a specific problematic area.

UV Irradiation

Ultraviolet (UV) irradiation is the preferred method for disinfection. This radiation acts on the epidermis and deeper layers of the skin, where it contributes to the production of protective agents, hormones and vitamins, in addition to its ability to kill bacteria.