CE Certification product



Aesthetic Technology Teeth Whitening Therapy

The best aesthetic beauty solutions right here!

✓ Dental enhancement

✓ Shade lifting

✓ Whitening

✓ Brightening



Adjustable LED triple color light with a wavelength of 460nm-490nm. The three colors represent different effects. Red light is for desensitization. Blue light is for whitening. Purple light is for sterilization

Blue, Red or purple light density of 200'300mw/cm2 assists with teeth discoloring due to smoking, aging, food, and health habits.

Each session lasts 10 – 15 minutes, but machine can operate 30 minutes at maximum setting.

According to dental sensitivity, choose the power strength between 40% and 100% for comfort.

This cold light operation system comes complete with two sized reusable mouth pieces, and Protective glasses

System
Specifications:

| Voltage | 110/220V and 60Watt(Max) |
|--------------------------|---|
| Dimensions | Length 130cm Machine, 130cm arm at most. Feet width 59 cm |
| Treatment time | 10 - 15 Minutes |
| Blue light wavelength | 460- 490 nm |
| Blue light density | 200- 300MW/cm2 |
| Weight | 18kg |







